**Ing. Josef Vašíček**

tel.: +420 723 169 089

Hrádek 60, 530 45 Opatovice n. Labem

e-mail: [vasicek.statika@seznam.cz](mailto:vasicek.statika@seznam.cz)

Objednatel : **Chrudimská nemocnice, a.s.**  
Václavská 570, 537 27 Chrudim  
IČ: 27520561  
DIČ: CZ27520561

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
## 3-2. System Layout

---

To maximize the performance obtained from equipment included in the system and increase the overall performance of the system, not only are the equipment directly introduced to the examination room and control room important, but the auxiliary and relevant peripheral equipment in both rooms are important as well.

These must be completed before beginning installation of the equipment delivered from Shimadzu. Therefore, after carefully considering the layout and peripheral equipment requirements with the customer, work directly and closely with the system engineer at the customer's contractor to make sure all the processes and construction work is completed necessary to ensure the necessary auxiliary work.

This system requires a minimum of about 30 m<sup>2</sup>, divided into two separate rooms. Divide the space into an examination room and control room, prepared with a utility pit (200 mm wide by 150 mm deep) for running cable wiring to each room. If a pit cannot be provided, then cables must be routed along the floor, which requires providing steel cable covers. The next section describes precautions for laying out equipment in each room.

---

## 3-3. Examination Room

---

This is the room where patient imaging occurs. It requires plenty of space for patients to enter and exit the room. It also requires storage space for contrast media and other chemicals or apparatus, and a place for dressing.

---

### 3-3-1. Floor of Examination Room

---

Floor Strength:

The heaviest item in the FLEXAVISION system is the X-ray table (ZS-5D/DS), which is constantly subject to vibration and other movement from operating the equipment. Therefore, floor strength is an important consideration.

1. Always anchor the system to the floor with anchor bolts.
  - If the unit is not secured to floor or the floor is not sufficiently strong, vibration will remain during imaging and prevent obtaining quality X-ray images. Furthermore, it could also cause damage or major injury if the unit tipped over in an earthquake or other natural disaster.
2. If examination room floor strength requirements exceed the maximum weight limit per area for the floor, the floor must be reinforced. In particular, avoid the use of wooden floors.
  - If installing the system on a wooden floor, reinforce the floor as necessary.

(ZS-5D) Dead Weight  $\approx 9.8 \times 76.5 \text{ N (750 kgf)}$

Dead Weight =  $1.25 \times \text{Operating/Service Load}$

Operating/Service Load  $\approx 600 \text{ kg}$

Foot Print Area of System =  $0.63 \text{ m}^2$

Unit Load =  $0.75 \text{ t} / 0.63 \text{ m}^2 \quad 1.19 \text{ t} / \text{m}^2 \approx 0.12 \text{ kg} / \text{cm}^2$

(ZS-5DS) Dead Weight  $\approx 9.8 \times 93.1 \text{ N (912.5 kgf)}$


Dead Weight =  $1.25 \times \text{Operating/Service Load}$

Operating/Service Load  $\approx 730 \text{ kg}$

Foot Print Area of System =  $0.40 \text{ m}^2$

Unit Load =  $0.9125 \text{ t} / 0.40 \text{ m}^2 \quad 2.28 \text{ t} / \text{m}^2 \approx 0.23 \text{ kg} / \text{cm}^2$

The dead load factors in instrument vibration.

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2,0 x 2,1 . 2

8

If floor strength is below this value, the floor must be reinforced.

---

#### **Strength of Concrete**

Typically the compressive strengths of floors are between 100 kg/cm<sup>2</sup> and 400 kg/cm<sup>2</sup>, but more commonly between 150 kg/cm<sup>2</sup> and 250 kg/cm<sup>2</sup>. Tensile strength is about 10% and bending strength about 15% to 20% of compressive strength.

(JSME Mechanical Engineers' Handbook 5th Edition)

---

### **3-3-2. Ceiling of Examination Room**

---

To ensure maximum performance from the R/F system, a minimum ceiling height of 2550 mm is required. However, limiting the movement range of the imaging system allows installing the system in a room with a ceiling height of 2400 mm, though this prevents maximizing performance. (For an explanation of how to limit the movement range, refer to 3.1.2 User Configuration in Part 2 of the ZS-5D/DS Instruction Manual.)

---

### **3-3-3. Entrance/Exit of Examination Room and Equipment Access**

---

An entrance/exit opening of about 1800 mm wide by 2200 mm high is required not only for patient access, but also for gurneys, wheelchairs, and other equipment. If the opening will be used to move instruments to the site, provide at least the minimum height and width indicated below.

110 cm (width) x 180 cm (height) min. (ZS-5D)

130 cm (width) x 180 cm (height) min. (ZS-5DS)

If using an elevator, provide at least the minimum dimensions indicated below.

110 cm (width) x 220 cm (depth) x 220 cm (height) min. (ZS-5D)

130 cm (width) x 220 cm (depth) x 220 cm (height) min. (ZS-5DS)

---

### **3-3-4. Narrow Access to Examination Room (ZS-5D)**

---

For model ZS-5D, if the elevator or room access is too narrow to allow transporting the instruments as packed for shipping, equipment can be disassembled in two stages to the following dimensions, as shown in Figure 3-3.

- (1) For access widths 960 mm to 1100 mm,  
disassemble equipment to dimensions H930 mm x W950 mm.
- (2) For access widths 885 mm to 960 mm,  
disassemble equipment to dimensions H875 mm x W885 mm.